

Product name		0.1ml PCR Tube	0.2ml PCR Tube	0.1ml PCR Eight Tube	0.1ml PCR Eight Tube	0.1ml PCR Eight Tube	0.1ml PCR Eight Tube	0.2ml PCR Eight Tube	0.2ml PCR Eight Tube	0.2ml PCR Eight Tube	0.2ml PCR Eight Tube	0.1mL 96 well PCR plate(Transparent,No Skirt)	0.1mL 96 well PCR plate (Transparent,No Skirt,Print)	0.1mL 96 well PCR plate (Transparent,Full Skirt)	0.1mL 96 well PCR plate (Transparent,Full Skirt)	0.1mL 96 well PCR plate (Transparent,No Skirt)	0.1mL 96 well PCR plate (Transparent,No Skirt,Print)	0.1mL 96 well PCR plate (Transparent,Full Skirt)	0.1mL 96 well PCR plate (Transparent,Full Skirt)	0.1mL 96 well PCR plate (Transparent,No Skirt)	0.1mL 96 well PCR plate (Transparent,No Skirt,Print)	0.1mL 96 well PCR plate (White,No Skirt)	0.1mL 96 well PCR plate (White,Half Skirt)	0.1mL 96 well PCR plate (White,Full Skirt)	0.1mL 96 well PCR plate (Transparent frame, transparent tube,full skirt)	0.1mL 96 well PCR plate-BL (White frame, transparent tube,full skirt)	0.1mL 96 well PCR plate (Transparent,Half Skirt,AB)	0.1mL 96 well PCR plate (Transparent, no skirt, tearable)	0.1mL 96 well PCR plate (White,No Skirt,Tearable)	0.2mL 96 well PCR plate (Transparent,No Skirt)	0.2mL 96 well PCR plate (White,No Skirt)	0.2mL 96 well PCR plate (Transparent,Half Skirt)	0.2mL 96 well PCR plate (Transparent,Half Skirt,Print)	0.2mL 96 well PCR plate (Transparent,Half Skirt,AB)	0.2mL 96 well PCR plate (Transparent,Half Skirt,Removable)	0.2mL 96 well PCR plate (Transparent,No Skirt,Tearable)	0.2mL 96 well PCR plate (White,No Skirt,Tearable)					
Axygen®	MaxyGene II		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
Applied Biosystems®	2720		✓					✓	✓	✓	✓																				✓	✓	✓	✓	✓		✓	✓	✓			
	9600		✓					✓	✓	✓	✓																					✓	✓	✓	✓	✓		✓	✓	✓		
	9700		✓					✓	✓	✓	✓																					✓	✓	✓	✓	✓		✓	✓	✓		
	9800 Fast			✓	✓	✓	✓							✓	✓	✓	✓	✓	✓	✓	✓	✓						✓	✓	✓							✓	✓	✓	✓		
	Veriti®/Fast		✓		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
	7000		✓					✓	✓	✓	✓	✓																				✓	✓	✓	✓	✓		✓	✓	✓		
	7300		✓					✓	✓	✓	✓	✓																				✓	✓	✓	✓	✓		✓	✓	✓		
	7500/Fast	✓		✓				✓	✓	✓	✓	✓																	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	
	7700		✓					✓	✓	✓	✓	✓																		✓	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓
	7900HT/Fast			✓	✓	✓	✓								✓	✓	✓	✓	✓	✓	✓	✓						✓	✓	✓							✓	✓	✓	✓		
	StepOne™	✓		✓	✓	✓	✓																					✓	✓	✓									✓	✓	✓	
	StepOnePlus™	✓		✓	✓	✓	✓								✓	✓	✓	✓	✓	✓	✓	✓						✓	✓	✓								✓	✓	✓		
	ProFlexs™PCR System		✓					✓	✓	✓	✓	✓																					✓	✓	✓	✓	✓		✓	✓	✓	
	QuantStudio™ 12K	✓	✓					✓	✓	✓	✓	✓																✓	✓					✓	✓	✓	✓	✓		✓	✓	✓
	ViiA™ 7/Fast	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	3100																																✓	✓	✓	✓	✓	✓	✓	✓	✓	
	3130																																	✓	✓	✓	✓	✓	✓	✓	✓	✓
	3700																																	✓	✓	✓	✓	✓	✓	✓	✓	✓
	3730/3730x																																	✓	✓	✓	✓	✓	✓	✓	✓	✓
Eppendorf®	Mastercycler		✓					✓	✓	✓	✓									✓	✓	✓	✓	✓	✓												✓	✓	✓	✓	✓	
	Mastercycler ep		✓					✓	✓	✓	✓															✓	✓	✓										✓	✓	✓	✓	✓
	Mastercycler M384																																									
	Mastercycler nexus		✓					✓	✓	✓	✓																											✓	✓	✓	✓	
	Mastercycler nexus eco		✓					✓	✓	✓	✓																											✓	✓	✓	✓	
	Mastercycler nexus flat		✓																																							
	Mastercycler nexus flat eco		✓					✓	✓	✓	✓																															
	Mastercycler nexus gradient		✓					✓	✓	✓	✓																											✓	✓	✓	✓	
	Mastercycler nexus gradient eco		✓					✓	✓	✓	✓																											✓	✓	✓	✓	
	Mastercycler nexus CSX1		✓						✓	✓	✓																												✓	✓	✓	✓
	Mastercycler nexus CSX1e		✓						✓	✓	✓																												✓	✓	✓	✓
	Mastercycler nexus SX1		✓						✓	✓	✓																												✓	✓	✓	✓
	Mastercycler nexus SX1e		✓						✓	✓	✓																												✓	✓	✓	✓
	Mastercycler pro 384																																									